



Material			
No.	Plan	Article	Quantity
1	YP05103	2899400	2
2	XP01354	2992393	4
3	XP01346	2992391	2
4	XP01353	2992927	2
5	XP01350	2899900	2
6	D4x10 DIN 1481	2992925	2
7*	---	US00011	1

| \* - Not depicted




A (1:4)



A (1:4)  
(After assembly)

1911982 - complete

EP00200 / I

<b>SIDEBARS 969</b>				Plan No. EP00200 / I
(Luna 2   Luna Basic 2   Club   Club Vario)				
Scale: 1:8	Drafted	Verified	Confirmed	
Date	25/03/2015	26/06/2024	26/06/2024	
Initials	AZ	PM	TK	