

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

HORIZONTAL TRANSPOT SYSTEM / CLUB 1/2			
Scale: 1:13	Drafted	Verified	Confirmed
Date	2023-02-15	2023-08-28	2023-08-28
Initials	AZ	OY	AZ

Plan No.
E00052584/B.1



No.	Plan	Article	Qty 1 set	Qty 2 sets
1	P00052570DP	2672737	2	4
2	P00052573	2672637	0	2
3	YP05175	2704437	2	4
4	Wheel D75 No Brake	1001255	2	4
5	M10x55 DIN933	1005636	4	8
6	Wheel D75 Brake	1995807	2	4
7	M10 DIN985	1001683	4	8
8	Screw 5.3x25	H016007	4	8
9*/#	P00066288	US00022		
9*/##	P00066291	US00023	1	2
9*/###	P00066292	US00024		
10	---	---	-	-
11	---	---	-	-
12*	Cardboard top spacer	2998024	1	2
13**	XP00116	2990074	4	8

* - Not depicted

** - Use Sterknops from foot side frame / only for Club

- Club with wooden panels

- Club without wooden panels, Club Vario with wooden panels, Club UL


- Club Vario without wooden Panels

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

HORIZONTAL TRANSPOT SYSTEM / CLUB ^{2/2}

Scale:	1:13	Drafted	Verified	Confirmed
Date	2023-02-15	2023-08-28	2023-08-28	
Initials	AZ	OY	AZ	

Plan No.
E00052584/B.1



Refer to SOP S-PL.80.09 for packing instructions.