



Material IN			
No.	Plan	Article	Qty
1	P00030194	1100280	1
2	E00031155	---	1
3	E00031309	---	2
4	139I Pen Click	1918595	2
5	Plug in Tube D22	1990648	4
6	M6x30 ISO 7380	1009588	2
7	X000862	1006463	6
8	---	1001554	4
9	X001016	1001682	2
10	P00031518	3995009	2
11	139I Seat Cushion	1919029	1
12	139I Toilet Seat	1919028	1
13	139I COMMODE BUCKET	1918593	1
14	139I Left Footrest	1919030	1
15	139I Right Footrest	1919031	1
16	139I Complete Footrest	1919032	1
17	Star Knob	2990024	2
18	M6x30 ISO 7380	1009588	4
19	M6DIN125	1001564	4
20	M6DIN985	1001680	4
21	---	1918642	2
22	---	1918643	2
23	P00031516	3995013	1
24	X001994	1007053	2
25	P00036097	3995011	1
26	P00031370	3995001	8

D1.1	Selvam E	E00031309 was modified	26-07-2024	AR
E1.1	David Van Branteghem	1918642 was 1903448; 1918643 was 1903450	20/09/2023	TT
C1.1	David Van Branteghem	1001554 was 3900974	20/09/2023	SD
B1.1	David Van Branteghem	BOM modified	22/02/2023	SS
No.	Requested by	Change Description	Date	Init.

1918599			
139I NIGHT CHAIR			
Scale: 1:10	Drafted	Verified	Confirmed
Date	30-04-2022	26-07-2024	26-07-2024
Initials	AR	MP	DB
			Plan No. E00031635/E1.1 